AeroEclipse[®] II


BAN[®] Nebulizer

Trudell Medical International
725 Baransway Drive, London, Ontario, Canada N5V 5G4
+1-519-455-4862
info@trudellmed.com www.trudellmed.com

REF 105501

105827-001 C. 2024-06-12. *trade-marks and registered trade-marks of Trudell Medical International (TMI).
† trade marks and registered trade marks of their respective companies. © TMI 2022, 2024. All rights reserved.

INTENDED USE

The nebulizer is intended to be used to deliver aerosol drug to the patient's lungs. It is to be used by infant, child and adult patients. The patient is to be awake and can breathe unassisted. The environments of use include hospitals and clinics. The nebulizer is single patient, reusable.

Breath actuation only produces aerosol when the patient inhales. Aerosol is not produced when the patient exhales or takes a break in treatment. The clinical benefit is more drug is available to the lungs. Less drug is lost to the environment.

⚠ Caution: Breath actuation delivers drug on inspiration only. This may result in an increase of the dosage to the lungs. Health care providers should take this into consideration with regards to dosage.

If you have concerns regarding the use of the nebulizer, please speak with your health care provider.

Notes

- The nebulizer must only be used by a person who understands the instructions for use and can use the device safely. The treatment must be carried out by a responsible person if the patient is not able to use the nebulizer on their own. Infants, children, and persons with limited capabilities must be supervised.
- The nebulizer is designed to aerosolize solutions and suspensions.
- Only use with drugs prescribed by a health care provider.
- Do not use in an anaesthetic or a ventilator breathing system.

BEFORE EACH USE

Ensure these instructions have been read and are kept available at all times. Carefully examine the device. If there is any visible dust or other debris inside, then clean it. Check inside the **EZ Twist[®]** Tubing (**F**) end for fraying or debris, clean it out if needed. Replace the nebulizer and the tubing (**F**) immediately if they are damaged or have missing parts.

⚠ Warning: Do not leave the nebulizer unattended with children. Contains small parts that can cause choking and tubing that can cause strangulation.

INSTRUCTIONS FOR USE

1 Ensure the exhalation valve (**B**) on the mouthpiece (**A**) faces down.

⚠ Warning: It is important to prevent exhaled breath from blowing near the eyes. Refer to the drug product monograph for more information.

ⓘ Note: For patients unable to use the mouthpiece, the nebulizer can be used with a **ComfortSeal[®]** Mask.

2 Inhale and exhale through the nebulizer mouthpiece (**A**) to verify the movement of the green feedback button (**C**) on the top of the nebulizer. The button should move up and down freely.

Notes:

- To operate the nebulizer in continuous mode (**I**), rotate the mode selector (**G**) clockwise to the continuous mode position .
- To return the nebulizer to breath actuated mode (**H**), rotate the mode selector (**G**) counterclockwise to the breath actuated position .

3 Unscrew and remove the top of the nebulizer. Place prescribed drug into the nebulizer cup (**E**). Reattach the top and gently hand tighten.

ⓘ Note: Keep the nebulizer in an upright position so the drug does not spill.

4 Attach one end of the tubing (**F**) to the fitting in the bottom of the nebulizer and the other end to the flowmeter or compressed air source. Ensure both ends are securely attached.

5 Turn on the compressed air source. Set the flowmeter to 7 to 8 L/min with a compressed air source capable of delivering 345 kPa (50 psi).

6 Place the mouthpiece (**A**) in mouth, and inhale slowly and deeply. As you inhale, the green feedback button (**C**) will move into the down position. This indicates the nebulizer is producing aerosol. Exhale normally. As you exhale, the green feedback button (**C**) will move into the up position. This indicates the nebulizer is not producing aerosol.

Notes:

- Do not place the lips over the exhalation valve (**B**) on the bottom of the mouthpiece (**A**). Doing so will prevent the valve from functioning properly.
- Do not cover the inhalation ports (**J**) on the nebulizer. Doing so will reduce volume of inhaled air.
- Ensure the nebulizer remains upright for proper nebulization.

7 Continue to inhale from the nebulizer until it starts to sputter, or as instructed by your health care provider.

8 Turn off the compressed air source.

ⓘ Note: Disconnect the tubing (**F**) from the compressed air supply after use.

CLEANING

⚠ Caution: Follow the cleaning instructions. Failure to do so may affect drug delivery.

Nebulizer

The nebulizer (excluding the tubing (**F**)) can be cleaned after each use.

- The mode selector (**G**) should be in the breath actuated position .
- Remove the tubing (**F**) from the bottom of the nebulizer.
- Unscrew and remove the nebulizer top (**D**).
- Shake out excess drug.
- Use any of the following methods to clean the nebulizer.

A Hand Washing

- Soak the three parts in a mild solution of liquid dish detergent and lukewarm clean water for about 5 minutes. Agitate gently.
- Rinse thoroughly with clean water.
- Shake out excess water.

B Sterile Water

- Thoroughly rinse the three parts in sterile water.
 - Shake out excess water.
- 6** Briefly connect the tubing (**F**) to the bottom of the nebulizer cup (**E**) and the compressed air supply to clear the air pathway.

7 Allow the parts to air dry thoroughly before reassembling.

8 Reattach the nebulizer top (**D**), making sure to line up the white cylinder with the stem of the nebulizer cup (**E**). Gently hand tighten until a stop is detected.

9 Reattach the mouthpiece (**A**) to the nebulizer with the exhalation valve (**B**) facing down.

10 After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

EZ Twist[®] Tubing

To remove condensation on the inside of the tubing (**F**), connect one end of the tubing (**F**) to the compressed air supply. Allow the air to run without the nebulizer attached. This will force air through the tubing, drying the inside. To clean the exterior, wipe with a damp cloth. Allow to air dry thoroughly.

DISINFECTION

The nebulizer (excluding the tubing (**F**)) can be disinfected each day.

Follow the cleaning instructions. Do not dry or reassemble the nebulizer.

1 Use any of the following methods to disinfect the nebulizer.

A Isopropyl Alcohol

- Soak the three parts in 70 % isopropyl alcohol for 5 minutes.
- Rinse the parts thoroughly with sterile water.
- Allow the parts to air dry thoroughly before reassembling.

B Cidex[®]

- Follow the disinfectant instructions for use.
 - Allow the parts to air dry thoroughly before reassembling.
- 2** After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

Notes

- The nebulizer (excluding tubing) can be cleaned and disinfected up to 28 times.
- No modifications of the nebulizer are permitted.
- The nebulizer parts comply with the relevant requirements of ISO 27427.
- The materials used for the components may not be compatible with emulsions. In such cases, information should be sought from the drug supplier.
- Not made or manufactured with bisphenol A (BPA), di(2-ethylhexyl) phthalate (DEHP), latex or lead.
- This nebulizer can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State, the MHRA or the Saudi Food and Drug Authority, as appropriate.

Type	Pneumatically powered jet nebulizer	
Breath Actuation Operation	Triggered mechanically (by the patient on inhalation) by an inhalation flow over 15 L/min	
Driving Gas	Air	
Driving Gas Flow (at nominal 345 kPa (50 psi) pressure when connected to the nebulizer)	Maximum/Nominal: 8 L/min Minimum: 7 L/min	
Maximum Pressure for Tubing at 22 °C ± 2 °C Ambient Conditions	345 kPa (50 psi)	
Maximum Fill Volume	6 mL	
Maximum A-Weighted Sound Pressure Level	63 dBA	
Mass of Nebulizer (without tubing)	0.03 kg (0.07 lb)	
Operating Temperature Limits	+15 °C (59 °F) to +40 °C (104 °F)	
Operating Relative Humidity Range	15 % to 95 % RH	
Storage Temperature Range	-40 °C (-40 °F) to +60 °C (140 °F)	
Storage Relative Humidity Range	Ambient to 95 % RH	
Recommended Replacement Period	7 days	

ISO 27427 PERFORMANCE DATA	At 7 L/min	At 8 L/min
Aerosol Output ¹	14.8 mg	16.1 mg
Aerosol Output Rate ¹	3.4 mg/min	3.9 mg/min
Percentage of Fill Volume Emitted Per Minute ¹	6.9 %	7.8 %
Residual Volume ¹	1 mL	1 mL
Mass Median Aerodynamic Diameter (MMAD) ²	3.4 µm	
Geometric Standard Deviation (GSD) ²	2.3	
Particle size % > 5 µm ²	33 %	
Particle size % 2 µm to 5 µm ²	40 %	
Particle size % < 2 µm ²	27 %	
Respirable Fraction ² (5 µm and smaller)	67 %	

¹ 2 mL fill volume of sodium fluoride 2.5% (M/V) solution. TMI data on file.

² 3 mL fill volume of sodium fluoride 2.5% (M/V) solution. TMI data on file.

Notes

- A suspension and/or high viscosity solution can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output and/or aerosol output rate which then can be different from those disclosed in these instructions.
- Disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for paediatric or infant populations.

P/N

Part number

Qty

Quantity

UDI

Unique device identifier

Nebulizer

Nebulizer

Single patient – multiple use

Designed and developed in Canada

For optimum performance
replace after 7 days

AU SPONSOR

Australian Sponsor

Scan or contact Trudell Medical International to obtain a glossary of the symbols used on packaging and instructions for use.

AeroEclipse^{MC} II

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Trudell Medical International
725 Baransway Drive, London, Ontario, Canada N5V 5G4
+1-519-455-4862
info@trudellmed.com www.trudellmed.com/ca/fr-ca

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