

- en **Instructions for use/Technical description**
Permanent YASARGIL aneurysm clips
- USA Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
- de **Gebrauchsanweisung/Technische Beschreibung**
Permanente YASARGIL Aneurysmen-Clips
- fr **Mode d'emploi/Description technique**
Clips d'anévrisme permanents YASARGIL
- es **Instrucciones de manejo/Descripción técnica**
Clips de aneurismas permanentes YASARGIL
- it **Istruzioni per l'uso/Descrizione tecnica**
Clip per aneurisma permanenti YASARGIL
- pt **Instruções de utilização/Descrição técnica**
Clipes permanentes de aneurisma YASARGIL
- nl **Gebruiksaanwijzing/Technische beschrijving**
Permanente YASARGIL aneurysmaclips
- da **Brugsanvisning/Teknisk beskrivelse**
Permanente YASARGIL aneurysmeclips
- sv **Bruksanvisning/Teknisk beskrivning**
Permanenta YASARGIL aneurysmclips
- ru **Инструкция по применению/Техническое описание**
Несъемные аневризматические клипсы YASARGIL
- cs **Návod k použití/Technický popis**
Permanentní aneurysmatické svorky YASARGIL
- pl **Instrukcja użytkowania/Opis techniczny**
Stałe zaciski tętniakowe YASARGIL
- sk **Návod na použitie/Technický opis**
Permanentné svorky pre aneurizmy YASARGIL
- tr **Kullanım Kılavuzu/Teknik açıklama**
Kesintisiz YASARGIL anevrizma klipsi
- zh **使用说明 / 技术说明**
YASARGIL 永久性脑动脉瘤夹
- ko **사용 설명서 / 기술 설명**
영구 YASARGIL 동맥류 클립


1

2

3

4

5

6

Aesculap® Permanent YASARGIL aneurysm clips

Legend

- 1 YASARGIL Phynox aneurysm clip, example of a straight clip
 - 2 YASARGIL titanium aneurysm clip, example of a straight clip
 - 3 YASARGIL aneurysm clip, correctly positioned
 - 4 YASARGIL aneurysm clip, incorrectly positioned
 - 5 YASARGIL aneurysm clip, incorrectly positioned
 - 6 YASARGIL aneurysm clip, incorrectly positioned
- A Test point
B Force
C Length of the jaw part

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

General notes

- In the following text, "permanent Aesculap YASARGIL aneurysm clips" are referred to as "aneurysm clips" or "clips."
- General risk factors associated with surgical procedures are not described in this documentation.

Intended use

The permanent aneurysm clips are intended for permanent occlusion of cerebral aneurysms. These permanent aneurysm clips are intended for single use only.

Indications

Cerebral aneurysms

Absolute contraindications

Permanent aneurysm clips are contraindicated for all applications except the permanent occlusion of cerebral aneurysms.

- ▶ Do not use in case of known sensitivity to foreign implant materials.

Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation: Medical or surgical conditions (e.g. comorbidities) which could hinder the success of the operation.

Risks, adverse effects and interactions

The following possible side-effects and interactions are referred to in the framework of the legal informational obligation:

- Infections
- Sensitization to the clip materials
- Nerve lesions
- Vascular injuries with the possible consequence of hematoma or bleeding
- Organ and other tissue injuries
- Disruptions to wound healing as well as other types of surgical complications.

The application of aneurysm clips involves the following severe risks:

- Shifting or breakage of the aneurysm clip
- Scissoring of the jaw parts
- Rupture of the aneurysm due to perforation, incomplete fitting of the clip jaw on the neck of the aneurysm.
- Reduction in the width of the vessel due to the clip being positioned near the edge of the vessel in large aneurysms
- Cerebrovascular spasms and sudden death

Other adverse (side) effects of surgical wound infections, and general surgical complications.

Application

Documentation

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the clips and their dimensions
- Location of intraoperative landmarks

Each patient must be informed comprehensively about the properties of aneurysm clips, as well as about the surgical procedure.

Labels

Each pack of clips contains additional labels showing the article number and the individual serial number of the clip. To facilitate the subsequent radiological examination, the patient's file (for the hospital) and the patient ID (for the patient) should be marked with these labels.

Patient ID

Important information concerning the implanted aneurysm clip and the operation can be noted down on the patient ID. To facilitate the postoperative radiological examination, each patient must receive a patient ID.

The patient ID can be ordered separately from Aesculap.

Product

CAUTION

Damage to the clip due to incorrect handling, restriction of the functionality and changing of the closing force!

- ▶ Only remove clips from their sterile packaging or storage for application.
- ▶ Only remove and apply with Aesculap applying clips for aneurysm clips.
- ▶ The clip must never be manipulated with the fingers.
- ▶ Do not repeatedly open and close clip.
- ▶ Discard clips that appear changed or show signs of damage, (e.g. incorrect position of the jaws, bent parts or discolorations).

Note

Aesculap recommends supplying the permanent aneurysm clips prior to operation in their unopened sterile packaging; in this way, damage to the aneurysm clips is avoided, their functionality is ensured and their correct closing force is maintained.

- ▶ Remove clips from double sterile packaging or make ready in the storage system.
- ▶ Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Immediately discard any damaged products.
- ▶ Carefully grip using the Aesculap applying forceps for aneurysm clips and remove, see Fig. 3. In doing so, ensure that the clip is not gripped as shown in Fig. 4/5/6 and observe the following points:
 - Clips must be used only with the appropriate size of applying or removing forceps (mini, standard, or long).
 - Phynox clips must only be used with Aesculap applying forceps for Phynox aneurysm clips.
 - Titanium clips must only be used with Aesculap applying or removing forceps for titanium aneurysm clips.
 - Applying or removing forceps supplied by other manufacturers must not be used.
- ▶ Apply clip and ensure that occlusion of the aneurysm or vessel is achieved and that the aneurysm clip sits correctly on the aneurysm neck during and after the implantation and sits tight with the blood vessel.
- ▶ Check the position of the clip and correct if necessary.
- ▶ If necessary, insert additional clips.

Safety advice

Users

- Appropriate clinical training as well as a theoretical and practical command of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.
- Aesculap is not responsible for complications caused by:
 - Incorrect indication or implant selection
 - incorrect surgical technique
 - Incorrect combination of implant components, particularly combination with components of other manufacturers
 - Exceeding the limitations of the treatment method or non-observance of essential medical precautions
- The user has a duty to consult with the manufacturer if the preoperative situation regarding use of the implants is ambiguous.

Product

- ▶ Do not use damaged clips.
- ▶ To avoid damage to the aneurysm clips:
 - Always treat the aneurysm clips with appropriate care.
 - Never open aneurysm clips with your fingers.
 - Avoid any manual and/or mechanical manipulation (e.g. with fingers or instruments) of the aneurysm clips.
- ▶ To exclude the possibility of a galvanic reaction, implant permanent aneurysm clips of different materials and/or from different manufacturers such that they can never touch.

Patients having Titanium and Phynox aneurysm clips can undergo MRI examinations subject to the conditions listed below:

- Static magnetic field of up to 3 Tesla
- Highest spatial gradient of the magnetic field of 720 Gauss/cm

The clips can warm up by a maximum of 2.2 °C after 15 minutes of continuous scanning under the above conditions.

The quality of the MR images may be affected by artifacts in the immediate vicinity of the implanted clip.

Note

For more information on MR compatibility, see TA015333.

Note

Stronger magnetic fields or an extension of the imaging area can lead to a significant increase in the incidence of positioning errors and artifacts.

Aesculap cannot accept any responsibility for aneurysm clips that are handled inappropriately or not in accordance with the present instructions for use.

Sterility, cleaning and disinfection, storage

DANGER

Risk to patients!

- ▶ The clips are intended for single use only.
- ▶ Discard clips that have been in contact with blood, tissue, or body fluids.

- The clips are packed in individually labelled protective packaging.
- The clips are gamma-sterilized.
- ▶ Store clips in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not prepare clips that are past their expiry date or whose packaging is damaged, see Validated reprocessing procedure.
- ▶ Ensure that the clips are not damaged in any way.
- The product can be cleaned and sterilized up to 100 times if undamaged and used with proper care, and if it has not come into contact with blood, tissue or body fluids. Any additional processing falls under the responsibility of the user.
- The life of the product is limited by damage, handling, storage and transport of the product.
- A careful visual inspection before use is the best opportunity to recognize a product that is no longer functional.
- The biological compatibility of the product may be affected by the accumulation of cleanser residues. Monitoring this is the task of the user.

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessor is responsible for this. The specified chemistry was used for validation.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com. The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

Processability is endangered by contact with blood, tissue, or body fluids!

- ▶ Discard clips that have been in contact with blood, tissue, or body fluids.
- ▶ Handle the clips with new gloves only.
- ▶ Keep storage devices covered or closed.

Restriction of the clips' closing force due to incorrect handling!

- ▶ Do not use a cleaning brush to clean the clips.
- ▶ Do not open the clips.

Manual cleaning/disinfection

Danger to the patient!

- ▶ The product must only be cleaned mechanically!

Mechanical cleaning/disinfection with ultrasonic pre-cleaning

Note

It is essential that the cleaning and disinfection device used is officially recognized or certified as effective (e.g. fulfillment of EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Pre-cleaning with ultrasound

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Ultrasonic cleaning	RT (cold)	>15	1	D-W	Dr. Weigert neodisher® SeptoClean*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

Phase I

▶ Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.

Phase II

▶ Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Dr. Weigert neodisher® SeptoClean application solution 1 %*
III	Neutralization	>10/50	1	FD-W	Dr. Weigert neodisher® Z stock solution 0.1 %
IV	Intermediate rinse	>10/50	1	FD-W	-
V	Intermediate rinse	>10/50	1	FD-W	-
VI	Thermal disinfection	90/194	5	FD-W	-
VII	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine**

D-W: Drinking water

FD-W: Fully desalinated water

*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

** The following cleaning and disinfection devices were used to verify cleanliness: Miele 7836 CD

- ▶ Place the clips on the tray for aneurysm clips FT008R or FT009R (avoiding rinsing blind spots). When doing so, place a maximum of 5 clips in one recess.
- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection and checking

- ▶ Check each individual clip.
- ▶ Discard and do not use clips that have the following features:
 - Signs of damage
 - Incorrect jaw position
 - Bent components
 - Misalignment
 - Dirt

Color coding may fade in the decontamination process.

- ▶ Discard and do not use clips if the color coding is no longer clearly visible.

Packaging

- ▶ Sort the product in the tray for aneurysm clips FT008R or FT009R.
- ▶ Package trays appropriately for the sterilization process in a sterile packaging conforming to EN ISO 11607-1 (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Sterilization

- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to EN 285 and validated according to EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes for prion inactivation
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store the clips in their sterile packaging or in a tray for aneurysm clips FT008R or FT009R in a dust-free area protected from chemical fumes and extreme air or temperature changes.
- For further information on the appropriate storage trays recommended by Aesculap, please contact Aesculap or consult the brochures for titanium and Phynox aneurysm clips. These brochures can also be ordered from Aesculap.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Product description

The aneurysm clips are available in two different materials:

- Cobalt alloy (Phynox) according to ISO 5832-7
- Titanium alloy Ti6Al4V according to ISO 5832-3

The aneurysm clips are available in three different sizes (Mini, Standard and Long). Please contact Aesculap for more information and to order brochures on titanium and Phynox aneurysm clips.

The aneurysm clips are color-coded as follows in order to differentiate between size and application.

Designation	Size	Color-coded labelling
Phynox aneurysm clips	Mini Standard Long	Silver
Titan aneurysm clips	Mini	Jaw piece: Silver Spring: Purple
	Standard	Jaw piece: Silver Spring: Blue

Each aneurysm clip has an individual serial number.

Each pack of permanent aneurysm clips contains additional labels showing the article number and the individual serial number of the clip.

Each aneurysm clip has been packed and sterilized individually (dose min. 25 kGy).

Each aneurysm clip in its original packaging is delivered in a double sterile packaging, together with instructions for use and labels.

The closing force of each aneurysm clip is measured individually and is stated on the packaging. The closing force is measured at 1/3 of the length from the tip of the jaw part from the test point A in the middle of the contact surface. The Phynox aneurysm clips are measure with the jaw open to 0.5 mm, see Fig. 1. The titanium aneurysm clips are measure with the jaw open to 1 mm, see Fig. 2.

Note

Permanent clips that are suitable for implantation have a silver colored jaw area (titanium clips) or are colored completely silver (Phynox clips) in order to distinguish them from temporary clips that are not suitable for implantation.

The temporary aneurysm clips in the proximal spring area can also be marked with ring-shaped, black markings.

Applying and removing forceps for aneurysm clips are labeled by size (Mini, Standard or Long) and clip material (titanium, Phynox).

The applying and removing forceps for titanium aneurysm clips may also be labeled with the appropriate color for the corresponding size of clips.

For further information on appropriate applying and removal forceps, as well as for brochures for titanium and Phynox aneurysm clips, please contact Aesculap.

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