SuperNO₂VA™ Nasal PAP Ventilation System—Medium, Large


Instructions for Use

REF	Catalogue number	Rx ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Manufacturer
X	Quantity		Do not re-use
	Not made with natural rubber latex		Exp. Date
	Consult Instructions for Use	MR	MR Safe
	Caution		Not made with DEHP
MD	Medical Device		

en **English**

The SuperNO₂VA Nasal PAP Ventilation System is comprised of the following components:

1. Corrugated tubing with bag adapter
2. SuperNO₂VA Nasal PAP Ventilation Device
3. Elbow circuit
4. Strap
5. Oxygen tubing
6. Nasal mask O₂ cap
7. Manometer port cap
8. Manometer tubing
9. 2L hyperinflation bag

Manufactured by
Vyaire Medical, Inc.
 26125 North Riverwoods Blvd
 Mettawa, IL 60045 USA
 1-833-327-3284
 customersupport@vyaire.com
 www.vyaire.com

EC REP **EMERGO EUROPE** 2797
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands

© 2019 Vyaire. Vyaire, the Vyaire Logo and SuperNO₂VA are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates.

vyaire™

Intended Use

The SuperNO₂VA Nasal PAP Ventilation System has a mask that creates a seal when positioned over a patient's nose to direct oxygen and/ or air to the upper airway during the continuum of anesthesia care. To be used under clinical supervision with adequate alarms and safety systems for monitoring and treatment of ventilatory failure.

Indications for Use


The SuperNO₂VA Nasal PAP Ventilation System is intended for short-term (<24 hours) use on adult patients (>30 kg) during the continuum of anesthesia care. This is a single patient use disposable system.

Instructions for Use

Review the sizing guidance chart to assist in choosing the best fitting SuperNO₂VA Nasal PAP Ventilation Device. The clinician should confirm proper fit to achieve required seal prior to selecting.

L	Adult Large	
	Height Range	
	5' 10" - 7' 3"	
M	Adult Medium	
	Height Range	
	4' 7" - 6' 9"	

1. Remove SuperNO₂VA Nasal PAP Ventilation System from packaging.
2. Prior to the patient lying on the procedure table, place the head strap on the table head surface with the SuperNO₂VA Nasal PAP Ventilation Device off to the right side of the table.

3. Lay patient down on procedure surface with head on top of head strap.
4. Place the SuperNO₂VA Nasal PAP Ventilation Device on patient's face with the anesthesia circuit port pointing cephalad.

5. To secure the SuperNO₂VA Nasal PAP Ventilation Device to the patient's face, slide the left head strap into the head strap cleat and fold sideways over

the cleat. With downward force, apply the desired pressure to the SuperNO₂VA Nasal PAP Device onto the patient's face and secure the head strap. The head strap may be placed above, below, or directly over the ear, depending on the user.

6. Attach the pre-assembled oxygen tubing to an oxygen flowmeter or an oxygen source.
7. Uncap manometer port and connect and attach to appropriate manometer using tubing provided. Manometer port is located between adjustable flow control valve and the corrugated tubing.
8. Adjust gas flow to desired flow rate.
9. Turn adjustable flow control valve to the fully closed position.
10. Occlude patient connection knob on elbow by placing thumb over 15mm I.D. port opening and allow bag to inflate. Check for proper operation.
11. With elbow still occluded, adjust flow meter and flow control valve for desired flow and pressure.
12. Connect SuperNO₂VA Nasal PAP Ventilation Device to elbow of hyperinflation system. Expand corrugated tubing to desired length.

13. Secure all connections.
14. As bag inflates, initiate ventilation in accordance with recommended procedures. If bag overinflates, reduce gas flow and adjust flow control valve.
15. To provide supplemental oxygen: First, disconnect the anesthesia circuit from the circuit port. Second, remove the mask oxygen port cap and connect one side of oxygen tubing to the SuperNO₂VA Nasal PAP Ventilation Device's oxygen port and the other side to a transport oxygen tank or hospital wall supply oxygen.
16. To remove the SuperNO₂VA Nasal PAP Ventilation Device from the patient, simply detach the head strap from its connection and lift off the patient's face.

Contraindications

The SuperNO₂VA Nasal PAP Ventilation System is contraindicated for use in long-term ventilation conditions and treatment of sleep apnea.

Warnings

- Do not allow the SuperNO₂VA Nasal PAP Ventilation Device's seal to come into direct contact with the patient's eyes; such action could lead to an eye injury. Use of mask may result in: drying of the eyes, eye pain, eye infections, or blurred vision.
- Do not leave the head strap secured over the patient's ears for an extended period of time or use excessive pressure at any time; such action could lead to ear injury.
- The mask may not be suitable for those predisposed to aspiration.
- Discontinue using the mask if there is an allergic reaction to any part of the mask.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary depending on the pressure settings, patient's breathing pattern, mask selection, and the leak rate. Monitor periodically the patient's oxygen levels while system is in use.
- Prior to system use on a patient, clinician should ensure that all connections and valves are operating correctly by setting the desired flow and pressure of the system before connecting to the SuperNO₂VA Nasal PAP Ventilation Device. Failure to follow instructions could lead to patient injury.
- When the supplemental oxygen tubing is connected to the oxygen port and nasal ventilation is being performed, monitor peak airway pressures on the ventilator to avoid over-pressurization. Over-pressurization could cause lung injury (barotrauma).
- The SuperNO₂VA Nasal PAP Ventilation System shall be used in accordance with this Instructions for Use. Read all sections of this IFU prior to use. Improper use of this system may cause serious injury.
- Do not use the SuperNO₂VA Nasal Ventilation PAP System near spark or open flame.
- This is a single patient device and is only to be used by one patient. Do not clean or sterilize. Sterilization or use of cleaning solutions may result in retention of harmful residues or leave this device nonfunctional.
- Excess pressure may cause skin breakdown; evaluate the patient and periodically inspect for any signs of skin deterioration.
- Do not use the SuperNO₂VA Nasal PAP Ventilation Device with an oral airway; such action may lead to oxygen desaturation.
- This product should only be used in conjunction with a manometer to determine the appropriate delivery pressure.
- The SuperNO₂VA Nasal PAP Ventilation System is MRI safe when used with MRI compatible equipment.

Cautions

- The hyperinflation bag system is designed to be used by medical personnel familiar with manual ventilation and airway maintenance techniques.
- Do not autoclave, gas, or chemically sterilize the system.
- Re-use may degrade the performance of the product or contribute to cross contamination.
- Federal (USA) law restricts this device to sale by or on the order of a physician.

Technical Specifications

Deadspace Volume:

- Medium 63 mL
- Large 97 mL

Disposal Specifications

Dispose of the SuperNO₂VA Nasal PAP Ventilation System in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

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