
# AeroChamber Plus Flow-Vu<sup>®</sup> HOSPITAL

Anti-Static Valved Holding Chamber  
MOUTHPIECE / LARGE MASK

**INTENDED USE**  
This chamber is intended to be used along with a metered dose inhaler to deliver aerosol medication to the patient's lungs.

Manufactured by:  
**Trudell Medical International<sup>®</sup>**  
725 Third Street,  
London, Ontario, Canada N5V 5G4  
✉ customerservice@trudellmed.com  
☎ +1-519-455-7060  
www.trudellmed.com

Distributed by:  
**Trudell Medical Marketing Limited**  
758 Third Street,  
London, Ontario, Canada N5V 5J7  
✉ tmm1@tmm1.com  
www.tmm1.com

CANADIAN DESIGNED  
AND MANUFACTURED

## DIRECTIONS FOR USE

**AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup> HOSPITAL** Anti-Static Valved Holding Chambers can be used directly out-of-package. Before use, ensure these instructions and the instructions supplied with the inhaler have been read and are kept available at all times.

<p>Before use, carefully examine the chamber. Replace immediately if any defect is noticed.</p>	<p>Remove cap(s) from inhaler and chamber.</p>	<p>Shake the inhaler immediately before use as per the instructions supplied with it.</p>	<p>Insert the inhaler into the Backpiece of the chamber.</p>	<p>Put mouthpiece into the patient's mouth. Ask the patient to close their lips around the mouthpiece to ensure an effective seal. The <b>Flow-Vu<sup>®</sup></b> Indicator only moves if the patient has a good seal.</p> <p>Apply mask to face and ensure an effective seal. The <b>Flow-Vu<sup>®</sup></b> Indicator only moves if the patient has a good seal.</p>	<p>Press the inhaler at the beginning of a slow inhalation as the <b>Flow-Vu<sup>®</sup></b> Indicator moves toward the patient. The patient should breathe in slowly and deeply through the mouth until a full breath has been taken, holding their breath for 5 – 10 seconds, if possible. Otherwise, have them keep their lips tight on the mouthpiece breathing normally 2 – 3 times through the chamber after the inhaler is pressed. Ask the patient to SLOW DOWN if you hear the <b>FlowSignal<sup>®</sup></b> Whistle sound. It means they are inhaling too quickly. Administer one (1) puff at a time.</p>	<p>Press the inhaler at the beginning of a slow inhalation as the <b>Flow-Vu<sup>®</sup></b> Indicator moves toward the patient. Use the <b>Flow-Vu<sup>®</sup></b> Indicator to count patient breaths. Maintain seal for 5 – 6 breaths after inhaler is pressed. Ask the patient to SLOW DOWN if you hear the <b>FlowSignal<sup>®</sup></b> Whistle sound. It means they are inhaling too quickly. Administer one (1) puff at a time.</p>	<p>Follow instructions supplied with the inhaler on how long to wait before repeating steps 3 – 6.</p>

## REPROCESSING INSTRUCTIONS FOR MULTI-PATIENT USE

**AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup> HOSPITAL** Anti-Static Valved Holding Chambers should be reprocessed following the instructions provided below. Reprocess this device between patients. For continuous use (exceeding one week) reprocess weekly.

		<p><b>3 Refer to equipment manufacturer's instructions for pre-cleaning and pasteurization cycles.</b></p> <p>a. Secure components in mesh bag and then place in metal baskets. Each disassembled device should be placed in a separate bag.</p> <p>b. The devices are now ready for the reprocessing cycle.</p> <p>c. At the completion of the pasteurization cycle, remove parts, shake off excess water and place into a sterile dryer for a four hour cycle (minimum). Place the components in the dryer so that water would naturally drain out.</p> <p>d. Remove components from the dryer.</p> <p>Throughout the drying cycle, ensure that components from the different devices do not get interchanged.</p>			
<p>Remove the Backpiece. To detach the Frontpiece, twist chamber as shown. Remove mouthpiece cap (if applicable).</p>	<p>Inspect components for foreign debris and other gross contamination from normal use. If required, remove debris and / or contamination using an enzymatic cleaner prior to pasteurization.</p>	<p>The effectiveness of this procedure has been proven using an Olympic Medical Pasteurmatic<sup>®</sup> Compact device. Pasteur SOL III enzymatic detergent was used in the 11 minute wash phase (38°C). The pasteurization phase duration was 30 minutes at a temperature of 74°C.</p>	<p>To reassemble, fit the Frontpiece on the end of the chamber and twist firmly until securely locked into position. For mouthpiece models, the protective cap should always be placed on the mouthpiece when the chamber is not in use.</p>	<p>Center the Alignment Feature on the Backpiece with the <b>Flow-Vu<sup>®</sup></b> Indicator, as shown. Press firmly to attach the Backpiece.</p>	<p>Draw a ✓ in the next box available on the chamber body, using a permanent marker that does not come off during reprocessing, to indicate a completed pasteurization cycle. Discard the chamber after 5 cycles. Based on validation testing, the <b>AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup> HOSPITAL</b> Chamber can be reprocessed up to five times.</p>

## STERILIZATION

- This chamber can be sterilized using low-temperature systems such as Steris<sup>®</sup> Vaporized Hydrogen Peroxide Sterilizer or Amsco<sup>®</sup> V-Pro<sup>®</sup> 1, V-Pro<sup>®</sup> 1 Plus or V-Pro<sup>®</sup> maX sterilizers. Remove backpiece prior to sterilization.
- The chamber can be sterilized in STERRAD<sup>®</sup> Hydrogen Peroxide Gas Plasma Technology: STERRAD 100NX (Standard and EXPRESS Cycle), STERRAD NX (Standard Cycle), STERRAD 100S (Short cycle), 200 (Short cycle), STERRAD 50.

## INSPECTION FOLLOWING REPROCESSING OR STERILIZATION

- Ensure that device is completely dry.
- No warping or distortion of any components is evident, including the **Flow-Vu<sup>®</sup>** Inspiratory Flow Indicator.
- Device printing is still legible.
- If any components appear damaged, the device should be discarded.

## Notes

- Chamber should be discarded after 5 reprocessing cycles. Environmental conditions, storage and frequency of use can affect product lifespan.
- This product contains no latex.
- When not in use, store in suitable container.
- If you have questions about the performance or usability of this product, please contact your TMM1 Sales Representative.

## ⚠ Cautions

- Do not leave the chamber unattended with children.
- Do not disassemble the chamber beyond what is recommended in the Reprocessing Instructions or damage may result.

- The **AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup> HOSPITAL** Chamber is a medical device. To ensure proper performance, it should only be reprocessed according to these instructions.
- Ensure directions for use have been read prior to use and are kept available at all times.

## ⚠ Warning

- This product has been designed for multi-patient use when reprocessed according to instructions and should therefore only be used within the hospital. This product is not intended for resale.
- Product may be permanently damaged if exposed to temperatures above 80° C.

**AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup> HOSPITAL** Chambers comply with CSA Standard Z264.1-02., Spacers and Holding Chambers for use with Metered Dose Inhalers.

101881-001 Rev C. \* trade-marks and registered trade-marks of Trudell Medical International. † trade-marks and registered trade-marks of respective companies. © TMI 2012, 2013. All rights reserved.

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