



Vyntus™ WALK  
mobile six-minute walk test

# Vyntus™ WALK

the complete mobile solution is not only comfortable for the patient to wear, it's also easy for technicians to use, and most importantly, it is key in standardizing your laboratory testing procedure.

**Vyntus™ WALK** lets patients perform a six-minute walk test (6MWT) wearing wireless sensors that communicate data to a tablet PC.

- The workflow-driven application seamlessly integrates with our SentrySuite® software platform for central data management and reporting.



## Benefit from:

- **Broad clinical value:** generates a diagnostic and prognostic measure of functional capacity and physical fitness.
- **Mobile solution:** works with an tablet PC and interfaces with the patient worn secure bluetooth oximeter.
- **Further data:** lets you manually enter blood pressure, oxygen supply, oxygen supply type and rate of perceived exertion.
- **Quality checks:** helps check the credibility of your data.
- **Sensor-type flexibility:** can be used with soft, flex, ear and forehead SpO<sub>2</sub> sensor types.



# Go with the flow

The intuitive Vyntus™ WALK tablet APP guides you through the standardized workflow. Accurate test data keeps you focus on our patient.




Handy pop-ups help you guide the patient

View real-time trends that display the patient's heart rate, SpO<sub>2</sub> and other events


Quick event buttons help you and your patient get to the next step

# Smart steps


Follow your patient, acquire and track the results, then verify and report all via your wireless tablet PC.




Baseline



Baseline SpO<sub>2</sub> Trend



Pause



Recovery



Baseline-Recovery




Result

# Vyntus™ WALK – Standalone Mode of Operation

Tablet based patient demographics entry - Testing - Direct PDF result


## Two simple ways to print your report

- WIFI direct print
- Print via USB cable



## One simple way to store the report PDF


- Network share drive /  
shared folder



# Vyntus™ WALK – SentrySuite Connected Mode of Operation

## Your benefits:

- Trend your patient
- Combined reporting like e.g. Spirometry or Diffusion
- Treatment efficiency (Pre/Post)
- Connect with WIFI or USB Tethering
- High data security during transfer with end to end encryption




# All good things come in threes

The Vyntus™ WALK is compatible with three sensor technologies – choose the sensor that best fits you and your patients' needs.

## Forehead sensor –

a convenient application site on the forehead for exercise testing

- Can be used when patients need their hands to support walking.
- Cleanable sensor – disposable sensor holder.



## Ear clip sensor –

offers flexibility for longer-term SpO<sub>2</sub> monitoring when digits are not the preferred application site

- Reusable ear clip solution.
- Can be used when patients need their hands to support walking.



## Soft clip sensor –

offers ultimate flexibility for SpO<sub>2</sub> measurement

- Offers superior patient comfort by providing an optimal distribution of pressure.
- Three sizes (S, M, L) to fit a variety of patient hands from children to adults.



All sensors are connected to the strap



# Vyntus™ WALK – Technical Specifications

## Parameters

Measured, typed in and calculated parameters

6MWD [m], 6MWW [kg\*m], 6MWS [m/s], Duration of Test [mm:ss], Lap Distance [m], Lap Count, Distance Final Lap [m], Number of Pauses, Duration of Pauses [mm:ss], Dyspnea [RPE], Exertion [RPE], SpO<sub>2</sub> [%], Lowest SpO<sub>2</sub> during Test [%], Ti88 [mm:ss], HR [bpm], Highest HR [bpm], BP Systolic [mmHg], BP Diastolic [mmHg], Supplied O<sub>2</sub> Flow [L/min], Supplied O<sub>2</sub> Type, Walking Aid, FEV1%Pred [%], BODE Index, BMI

## Tablet

Minimum requirements	Operating system	Android 4.4.x or 5.x + Recommended Android 7.x
	Processor	Intel® Atom™ (1.6 GHz Dual-Core or higher) or Quad-Core ARM (1,3 GHz or higher)
	Screen size	7" tablet screen
	Screen resolution	1280 x 800 pixels
	Bluetooth®	2.0 or higher
	Wireless type	WiFi 802.11n

## Network / WLAN

Minimum requirements	Connection speed	150 Mbit/s
	Printer	Wireless supporting printer for stand-alone mode

## SpO<sub>2</sub> device

Type	Nonin	WristOx2® Model 3150 Bluetooth® classic and BLE
Ranges	Oxygen saturation range (SpO <sub>2</sub> )	0 % to 100 %
	Pulse rate range	18 to 321 beats per minute (BPM)
Display features	SpO <sub>2</sub>	3-digit numeric
	Pulse rate	3-digit numeric, Pulse strength bar graph: 9 levels
	Battery status	4 levels (full, half, low, critical)
	Indicators	Sensor Fault, Bluetooth connectivity, SmartPoint, Pulse Signal
Internal Power	Battery	Two 1.5 Volt AAA alkaline batteries
Weight / Dimensions	Weight (approx. with strap and sensor)	71 grams (2.5 oz.)
	Dimensions (without strap and sensor)	56 mm × 74 mm × 20 mm (W × H × D) (2.20" × 2.91" × 0.79")

## SpO<sub>2</sub> Sensors

SpO <sub>2</sub> Accuracy (Arms*) @ 70 – 100% Adult/Pediatric		
Soft sensor (small, medium, large)	No motion	± 2 digits
	Motion	± 3 digits
	Low perfusion	± 2 digits
Flex sensor adult, child, infant Reflectance sensor	No motion	± 3 digits
	Motion	N/A
Ear clip sensor	Low perfusion	± 2 digits
Pulse Rate Accuracy (Arms*) Adult/Pediatric		
All sensor models	No motion @ 18 – 300 BPM	± 3 digits
	Motion @ 20 – 250 BPM	± 3 digits
	Low perfusion @ 40 – 240 BPM	± 3 digits

## Conditions

Operating conditions	Temperature	0 °C to +35 °C (+32 °F to 95 °F)
	Relative humidity	15 to 95 % RH, non-condensing
Storage/Transportation conditions	Temperature	-20 °C to +50 °C (-4 °F to 122 °F)
	Relative humidity	15 to 95 % RH, non-condensing




## NOTES

\*  $\pm 1$  Arms represents approximately 68 % of measurements.

## GLOBAL HEADQUARTERS

Vyairē Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045  
USA

 Vyairē Medical GmbH  
Leibnizstrasse 7  
97204 Hoechberg  
Germany



[vyaire.com](http://vyaire.com)

For EU, US, Saudi Arabia and United Arab Emirates distribution only.

All trademarks are property of their respective owners.

© 2019 Vyairē Medical, Inc. or one of its affiliates. All rights reserved. Vyairē, the Vyairē Medical logo, SentrySuite and Vyntus are trademarks or registered trademarks of Vyairē Medical, Inc. or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling VYR-GLB-1900001 (2.0).