


Filter Kit for use with **AeroEclipse®**

BAN™ Nebulizers



- ✓ Filters fugitive emissions (aerosolized medication)¹, viruses² and bacteria² from the patient's exhaled breath
- ✓ Reduces exposure of healthcare providers to potentially harmful emissions

The filter kit is a single patient, single use accessory to the **AEROECLIPSE® BAN™** Nebulizer.

Filter Kit, shown with **AEROECLIPSE® II BAN™** Nebulizer

Aerosol Performance:

Aerosol Output ¹ of AEROECLIPSE® BAN™ Nebulizer with Filter Kit compared to the BAN™ Nebulizer alone	≥ 90 %
Aerosolized Medication Emissions ¹ when assembled with AEROECLIPSE® BAN™ Nebulizer	≤ 1 %

Filtration Performance:

Viral Filtration Efficiency ²	≥ 99.99 %
Bacterial Filtration Efficiency ²	≥ 99.999 %

1. TMI, Data on file 30 August 2020. 2. Breathing Filter Product Report, 24 October 2020, on file TMI. 115010-001 Rev.C. ©/TM trade-marks and registered trade-marks of Trudell Medical International (TMI). © TMI 2019, 2020, 2021. All rights reserved.



aeroeclipse.com

