

# Aesculap<sup>®</sup> AdTec<sup>®</sup> single use

Single use bipolar instruments



Aesculap Endoscopic Technology






# Single use bipolar instruments

## Bipolare Einmalinstrumente



- Excellent coagulation results
- Scissors with outstanding cold cutting performance
- Jaw parts feature ceramic insulation
- Insulated to the tip
- Non-protruding joints
- Ergonomic handle design
- One-to-one tactile feedback
- Compatible with all electrosurgical units
- Shaft 360° rotatable
- Ausgezeichnete Koagulationsergebnisse
- Schere schneidet hervorragend auch ohne Strom
- Keramikisolation im Maulbereich
- Isolation bis an die Maulteile
- Keine austretenden Gelenkteile
- Ergonomische Griffform
- Hervorragende Taktilität
- Kompatibel zu allen Elektrochirurgie-Generatoren
- Schaft 360° drehbar



			Pack (pieces) Packung (Stück)	Diameter Durchmesser	Length Länge
PL700	 1/1	METZENBAUM single use scissors, tips curved to the left METZENBAUM Einmalschere, Blätter nach links gebogen	6	5 mm	310 mm
PL701	 1/1	MARYLAND single use grasping forceps, fenestrated MARYLAND Einmalfasszange, gefenstert	6	5 mm	310 mm
PL702SU	 1/1	MARYLAND single use grasping forceps MARYLAND Einmalfasszange	6	5 mm	310 mm
PL703SU	 1/1	Macro single use grasping forceps Macro Einmalfasszange	6	5 mm	310 mm
		<b>Bipolar single use cable</b> , suitable for <b>Bipolares Einmalkabel</b> , geeignet für		pin distance Pinabstand	Length Länge
GK290SU		Dolly, EMC, ERBE Int.		22 mm	4 m
GK291SU		Aesculap, BOWA, Lamidey, Söring, Valleylab		28.6 mm	4 m
GK292SU		Codman (CMC II, CMC III)		38 mm	4 m

**NEW**

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

The main product trademark 'Aesculap' and the product trademark 'AdTec' are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C09411

1111/2.0/5