



AeroChamber® mini Anti-Static Holding Chamber

AeroChamber® mini Anti-Static Holding Chamber is designed for the delivery of metered dose inhaler (MDI) medications in mechanical ventilation circuits or with a manual resuscitation bag. Compact, the chamber minimizes dead space.

Refer to the [instructions for use](#) for additional information.



- Small 110mL volume chamber, suitable for application in neonates and other patients with lower tidal volumes¹
- Anti-static improves the reliability of medication delivery²
- Ready to use directly out of the package – time saving with no pretreatment needed

Part Number	Case Quantity
875110010	10

SINGLE-PATIENT – MULTIPLE USE | Replace with the ventilator circuit. Do not clean or disinfect.

PACKAGING: Flow-wrapped with instructions for use.

Not made or manufactured with bisphenol A (BPA), phthalates, latex or lead.

AeroChamber®
mini

Anti-Static Holding Chamber

1. Nagel M, et al. *European Respiratory Journal* 2009;34(S53):P2036.

2. Suggett J, et al. *Journal of Aerosol Medicine and Pulmonary Drug Delivery* 2015;28(4):A-1-A-25.

Manufactured by:

+1-866-510-0004
info@trudellmed.com

trudellmed.com



Designed and Developed in Canada
Manufactured in Canada with USA and Canadian parts.

Canadian Hospital Distribution:
Trudell Healthcare Solutions
customerservice@trudellhs.com
www.trudellhs.com

