

vyaire™
MEDICAL



AirLife™
Noninvasive ventilation masks

AirLife™ noninvasive ventilation masks


Our noninvasive ventilation (NIV) full face and nasal masks offer clinicians a reliable option while providing patients comfort and security.

The adjustable four-point headgear enables ease in fit and adjustment while the flexible swivel elbow reduces torque applied to the seal between the mask and patient, helping to minimize pressure points.


Our masks feature a lightweight silicone design that contours to the patient's face for a comfortable fit, and the quick-release clip eliminates the need to refit after removal.

The odorless, easy-to-clean masks include both vented and non-vented options and are available in three sizes to fit a wide range of patients.

The mask flexes when the circuit is pulled.



Mask when flexed



Mask not flexed



Full face mask

Key features

- Flexible swivel connector enables ease in positioning of the circuit
- Comfortable four-point headgear with crown strap for better stability and sealing
- Adjustable forehead rest for comfort and stability
- Soft seal that conforms to the face, reducing need to overtighten
- Quick-release strap
- Made without DEHP

AirLife NIV masks (10/cs.)

| Cat no. | Product description |
|---------|--|
| NIV040S | Full face mask, vented with anti-asphyxia valve (AAV), small |
| NIV040M | Full face mask, vented with AAV, medium |
| NIV040L | Full face mask, vented with AAV, large |
| NIV043S | Full face mask, non-vented with AAV, small |
| NIV043M | Full face mask, non-vented with AAV, medium |
| NIV043L | Full face mask, non-vented with AAV, large |
| NIV041S | Full face mask, non-vented, small |
| NIV041M | Full face mask, non-vented, medium |
| NIV041L | Full face mask, non-vented, large |
| NIV042S | Nasal mask, vented, small |
| NIV042M | Nasal mask, vented, medium |
| NIV042L | Nasal mask, vented, large |
| NIV044S | Nasal mask, non-vented with AAV, small |
| NIV044M | Nasal mask, non-vented with AAV, medium |
| NIV044L | Nasal mask, non-vented with AAV, large |

For more information, please call Customer Service at **800.323.9088**.



GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA

[vyaire.com](https://www.vyaire.com)

U.S. ONLY For U.S. distribution only.

© 2020 Vyaire. Vyaire, the Vyaire logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-US-2000139

vyaire[™]
M E D I C A L