



# Butterfly iQ+/iQ Accessories





# Butterfly Accessories



## Butterfly Hard Case

275 x 265 x 88 mm

1.53 kg

*Plastic: PC-ABS (poly carbonate/acrylonitrile butadiene styrene).*

*Handle: TPU. Inner tray: Silicone.*



## Butterfly Soft Case


210 x 178 x 38 mm

280 g

*Materials: 1050D Ballistic Nylon.*



# Butterfly Accessories



## Extra Butterfly Charger

122 x 21 x 62 mm

0.05 kg

*Materialien: Polycarbonate.*



## Butterfly Holster

61 x 50 x 63 mm

0.08 kg

*Materials: Impact modified polycarbonate  
and durable silicone.*



# Butterfly Accessories



**USB-C 2.0 Cable**  
(iQ+/iQ+ Vet only)

**1.5 m**

1500 x 24 x 18 mm

41 g

*Materials: Durable TPE.*



**Lightning® Cable**  
(iQ+/iQ+ Vet only)

**1.5 m**

1500 x 24 x 18 mm

41 g

*Materials: Durable TPE.*



**USB-C 2.0 Cable**  
(iQ+/iQ+ Vet only)

**2.5 m**

2500 x 24 x 18 mm

74 g

*Materials: Durable TPE..*



**Lightning® Cable**  
(iQ+/iQ+ Vet only)

**2.5 m**

2500 x 24 x 18 mm

74 g

*Materials: Durable TPE.*