

EasyGoVac Aspirator

Features

Lightweight at only 9.6 lbs

2 hour battery runtime

Low battery light

Dual pump protection systems with float shutoff and inline hydrophobic filter safeguards

Can be used with any standard 800cc canister

Easy-to-read gauge

Replaceable inline hydrophobic bacteria filter

Product Specifications

Vacuum range: 51-483 mmHg

Vacuum mode:
Continuous

Height: 8.8 in (22.4 cm)

Width: 6.7 in (17.0 cm)

Length: 8.9 in (22.5 cm)

Weight: 9.6 lbs (4.4 kg)

Suction canister: 800cc
(Disposable with float shut-off)

Warranty

2-year on device

6-month on battery

Ordering Information

EasyGoVac (PM65)

Includes EasyGoVac with 800cc canister, hospital-grade power supply and battery pack with charger, suction tubing, inline hydrophobic bacteria filter, elbow connector and nylon carrying case



PM65

The EasyGoVac battery-powered aspirator combines ease of use with powerful suction to create a portable unit that is ideal for hospitals, nursing homes, ambulances and homecare. EasyGoVac produces vacuum up to 483 mmHg and provides dependable continuous suction power to remove respiratory fluids. The high-performance battery delivers 2 hours of runtime, more than any aspirator of comparable size. The “Low Battery Light” lets users know when the unit needs to be recharged. EasyGoVac is easy to clean and the closed design guards against bacterial growth. Safety features and safeguards in the pump prevent overflow and possible contamination.

Aspirator Accessories

Canisters

502519-10

800cc canister with lid (case of 10) - PM60 and PM65

1906-12

1200cc canister with lid (case of 12) - PM61 and PM63

502687-12

2000cc canister with lid (case of 12) - PM61 and PM63

502540

Canister bracket - PM61 and PM63

Roll Stands

PM60-1

Roll stand for use with any Precision Medical aspirator or compressor

PM60-2

Roll stand with canister bracket for use with any Precision Medical aspirator

Power Supply and Cord

503974

Battery charger power supply - PM65

502954

Auto DC power cord - PM65

Carry Bag

502956

Carry bag - PM65

Tubing and Filters

1690

Inline hydrophobic bacteria filter - PM60 and PM63

1955

Patient suction tubing (6 feet)

502012

Inline hydrophobic bacteria filter assembly - PM60

502391

Small air filter - PM63

502471

Large air filter - PM63

502690

Inline hydrophobic bacteria filter & tube assembly - PM61 and PM63

502698

Plastic tubing replacement kit (2" and a 12" section) - PM60

503293

Inline hydrophobic bacteria filter assembly with tubing and elbow - PM65

506047

Inline hydrophobic bacteria filter - PM65

