

Precision Medical offers a variety of flowmeter solutions to meet the needs of your customer. Please reference the charts on page 2 for the full list of product numbers.

Most Common Flowmeters:

Product Code:	Description:
1MFA1003	Oxygen, Chrome, 0-15 l/min, DISS Female H.T.
8MFA1003	Oxygen, Compact, 0-15 l/min, DISS Female H.T.



1MFA1003



8MFA1003



DISS Female Hand Tight

Most Common Accessories:

DISS "Y" Assemblies with Check Valves

Product Code:
7211 - Oxygen
7212 - Air
7213 - Vacuum
7214 - Nitrous Oxide
7240 - Nitrogen



Disposable Tubing Nipples

Product Code:
0119 - Clear
0120 - Yellow
0121 - Black
0122 - White
0123 - Green



There are 4 important questions to ask your customer:

1. What type of gas is required? (Oxygen or Air)
2. How much flow is required? (Flow Range in l/min or cc)

O2 0-200 cc	O2 0-1 l/min	O2 0-3.5 l/min	O2 0-5 l/min	O2 0-15 l/min	Air 0-15 l/min	O2 2-26 l/min	O2 3-35 l/min	O2 0-70 l/min	Air 0-70 l/min
----------------	-----------------	-------------------	-----------------	------------------	-------------------	------------------	------------------	------------------	-------------------

3. What type of fitting is required? (DISS, Ohmeda, etc.)



DISS Female Hex Nut



DISS Female Hand Tight



DISS Male



Ohmeda



Chemetron



Oxequip




Puritan - Bennett



Schrader



Medstar












TRUDELL
HEALTHCARE SOLUTIONS

4. What model of flowmeter is required? (Compact or Chrome)

Compact Flowmeters

Lightweight, Compact, Accurate, Dependable

	O2 0-8 l/min	O2 0-15 l/min	Air 0-15 l/min
No Adaptor (1/8 NPTF)	8MFA4001	8MFA1001	8MFA2001
 DISS Female Hex Nut	8MFA4002	8MFA1002	8MFA2002
 DISS Female Hand Tight	8MFA4003	8MFA1003	8MFA2003
 DISS Male	8MFA4004	8MFA1004	8MFA2004
 Ohmeda	8MFA4005	8MFA1005	8MFA2005
 Chemetron	8MFA4006	8MFA1006	8MFA2006
 Oxequip	8MFA4007	8MFA1007	8MFA2007
 Puritan- Bennett	8MFA4008	8MFA1008	8MFA2008
 Schrader	8MFA4009	8MFA1009	8MFA2009
 Medstar	8MFA400M	8MFA100M	8MFA200M

Compact Models:



8MFA4001



8MFA1001



8MFA2001

Chrome Models:



6MFA1001












4MFA1001




1MFA3001

Chrome Flowmeters

Durable, Accurate, Precise, Dependable

	O2 0-200 cc	O2 0-1 l/min	O2 0-3.5 l/min	O2 0-5 l/min	O2 0-15 l/min	Air 0-15 l/min	O2 2-26 l/min	O2 3-35 l/min	O2 0-70 l/min	Air 0-70 l/min
No Adaptor (1/8 NPTF)	6MFA1001	4MFA1001	1MFA3001	1MFA3501	1MFA1001	1MFA2001	1MFA1101	1MFA1401	1MFA8001	1MFA9001
 DISS Female Hex Nut	6MFA1002	4MFA1002	1MFA3002	1MFA3502	1MFA1002	1MFA2002	1MFA1102	1MFA1402	1MFA8002	1MFA9002
 DISS Female Hand Tight	6MFA1003	4MFA1003	1MFA3003	1MFA3503	1MFA1003	1MFA2003	1MFA1103	1MFA1403	1MFA8003	1MFA9003
 DISS Male	6MFA1004	4MFA1004	1MFA3004	1MFA3504	1MFA1004	1MFA2004	1MFA1104	1MFA1404	1MFA8004	1MFA9004
 Ohmeda	6MFA1005	4MFA1005	1MFA3005	1MFA3505	1MFA1005	1MFA2005	1MFA1105	1MFA1405	1MFA8005	1MFA9005
 Chemetron	6MFA1006	4MFA1006	1MFA3006	1MFA3506	1MFA1006	1MFA2006	1MFA1106	1MFA1406	1MFA8006	1MFA9006
 Oxequip	6MFA1007	4MFA1007	1MFA3007	1MFA3507	1MFA1007	1MFA2007	1MFA1107	1MFA1407	1MFA8007	1MFA9007
 Puritan- Bennett	6MFA1008	4MFA1008	1MFA3008	1MFA3508	1MFA1008	1MFA2008	1MFA1108	1MFA1408	1MFA8008	1MFA9008
 Schrader	6MFA1009	4MFA1009	1MFA3009	1MFA3509	1MFA1009	1MFA2009	1MFA1109	1MFA1409	1MFA8009	1MFA9009
 Medstar	6MFA100M	4MFA100M	1MFA300M	1MFA350M	1MFA100M	1MFA200M	1MFA110M	1MFA140M	1MFA800M	1MFA900M



TRUPELL
HEALTHCARE SOLUTIONS