

HEARTSYNC
Defibrillation • Cardioversion • Pacing • ECG Monitoring



36 MONTH SHELF LIFE
LONGEST ON THE MARKET

Extended shelf life reduces waste and saves money



60 INCH LEAD WIRES
LONGEST ON THE MARKET

Increased flexibility with equipment placement & clinician mobility



ADVANCED HYDROGEL

Unique gel provides maximum adhesion



SMALLER PAD SIZE

Allows for easier placement without crowding



LARGER CONDUCTOR PLATE

Larger conductive area disperses energy more effectively

COLOR CODED PACKAGING

The correct product, for your device.

VISIBLE LINER

Prevents slip hazard to clinician.

INSTRUCTIONS FOR USE

User checklist, on the release liner.

COLOR CODED PADS

Anterior electrode is white foam and the posterior electrode is blue foam. This helps train staff to look for colors to aid in pad placement.

LEADS OUT OPTION

Allows pads to be pre-connected prior to opening pouch

Latex-free



PHYSIO CONTROL




PHILIPS MEDICAL



CARDIAC SCIENCE



ZOLL MEDICAL




TRUDELL
HEALTHCARE SOLUTIONS

Defibrillation Electrodes

Defib Pads for all departments of the hospital



Vermed's HeartSync Defibrillation electrodes are designed with your needs in mind. Featuring the longest leadwires, 36 month shelf life, and ease of placement, they will meet the needs of your clinical staff.



PHYSIO CONTROL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
C100AC-PHYSIO	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
T100AC-PHYSIO	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
T100LOAC-PHYSIO	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Out	
PEDIATRIC-PHYSIO	Lifepak 9P, 15, 10C, 20, 12	In	

PHILIPS MEDICAL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
C100AC-PHILIPS	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
T100AC-PHILIPS	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
T100LOAC-PHILIPS	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Out	
PEDIATRIC-PHILIPS	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	In	

CARDIAC SCIENCE

Product #	For Defibrillator Models	Leads	Pad & Leadwires
T100-CS	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	Can be used in place of Cardiac Science Part #9131-001


ZOLL MEDICAL

Product #	For Defibrillator Models	Leads	Pad & Leadwires
C100AC-ZOLL	M-Series, E-Series, R-Series, X-Series	In	
T100AC-ZOLL	M-Series, E-Series, R-Series, X-Series	In	
T100LOAC-ZOLL	M-Series, E-Series, R-Series, X-Series	Out	
PEDIATRIC-ZOLL	M-Series, E-Series, R-Series, X-Series	In	

For additional information, defibrillator compatibility and usage visit: vermed.com/defib


PHYSIO CONTROL

Connector Style




PHILIPS MEDICAL

Connector Style




CARDIAC SCIENCE

Connector Style



ZOLL MEDICAL

Connector Style



758 Baransway Dr, London, ON N5V 5J7

1-800-265-5494 | www.trudellhs.com