



Infant Flow™ LP nCPAP System
Providing nCPAP solutions for your
most precious patient

Infant Flow™ LP Headgear

Part #	Description	Measurement	UM
777040XS	Infant Flow LP Headgear, Extra Small	17-21 cm	10 / package
777040S	Infant Flow LP Headgear, Small	21-26 cm	10 / package
777040SM	Infant Flow LP Headgear, Small/Med	24-28 cm	10 / package
777040M	Infant Flow LP Headgear, Medium	26-32 cm	10 / package
777040L	Infant Flow LP Headgear, Large	32-37 cm	10 / package
777040XL	Infant Flow LP Headgear, Extra Large	37-42 cm	10 / package

Infant Flow™ LP Bonnets

Part #	Description	Measurement	UM
777010	Infant Flow LP Bonnet, Size 000	18-20 cm	10 / package
777012	Infant Flow LP Bonnet, Size 00	20-22 cm	10 / package
777014	Infant Flow LP Bonnet, Size 0	22-24 cm	10 / package
777016	Infant Flow LP Bonnet, Size 1	24-26 cm	10 / package
777018	Infant Flow LP Bonnet, Size 2	26-28 cm	10 / package
777020	Infant Flow LP Bonnet, Size 3	28-30 cm	10 / package
777022	Infant Flow LP Bonnet, Size 4	30-32 cm	10 / package
777024	Infant Flow LP Bonnet, Size 5	32-34 cm	10 / package
777026	Infant Flow LP Bonnet, Size 6	34-36 cm	10 / package
777028	Infant Flow LP Bonnet, Size 7	36-38 cm	10 / package

Infant Flow™ LP nCPAP Generator / Circuits

Part #	Description	UM
7772011	Infant Flow LP nCPAP Generator/Circuit Kit (F&P MR 730 & 850 Humidifiers)	10 / case
7772021	Infant Flow LP nCPAP Generator/Circuit Kit (F&P MR 730& 850 Humidifiers)	20 / case

Infant Flow™ nCPAP Starter Kit

Part #	Description	UM
7772010	Infant Flow LP nCPAP Starter Kit	Each

Infant Flow™ LP Nasal Prongs

Part #	Description	UM
777000XS	Nasal Prongs, Extra Small	10 / box
777000S	Nasal Prongs, Small	10 / box
777000M	Nasal Prongs, Medium	10 / box
777000L	Nasal Prongs, Large	10 / box
777000XL	Nasal Prongs, Extra Large	10 / box

Infant Flow™ LP Nasal Mask

Part #	Description	UM
777002XS	Nasal Masks, Extra Small	10 / box
777002S	Nasal Masks, Small	10 / box
777002M	Nasal Masks, Medium	10 / box
777002L	Nasal Masks, Large	10 / box
777002XL	Nasal Masks, Extra Large	10 / box

Infant Flow™ LP Nasal Mask

Designed for maximum therapeutic effectiveness and patient comfort. The flexible bellow allows the mask to find a natural position over the infant's nose, minimizing leaks to maintain a consistent nCPAP level. An assessment window enables visualization of the infant's nares and septum. The five masks are color coded for easy size identification. Individually packaged 10 per box.



Part #	Description
777002XS	Nasal Masks, Extra Small
777002S	Nasal Masks, Small
777002M	Nasal Masks, Medium
777002L	Nasal Masks, Large
777002XL	Nasal Masks, Extra Large

Infant Flow™ LP Nasal Prongs

Designed for maximum therapeutic effectiveness and patient comfort, the prongs have high-profile design to minimize skin contact and provide visual verification of proper placement. The gently flared prongs help to form an effective seal and minimize leaks. The five prongs are color coded for easy size identification. Individually packaged 10 nasal prongs per box.




Part #	Description
777000XS	Nasal Prongs, Extra Small
777000S	Nasal Prongs, Small
777000M	Nasal Prongs, Medium
777000L	Nasal Prongs, Large
777000XL	Nasal Prongs, Extra Large

Infant Flow™ LP nCPAP Generator/Circuit Kit

The new generator, with dual-jet variable flow technology, effectively reduces the imposed WOB to maintain a constant CPAP level. The lower driving pressure reduces pressure back up into the humidifiers auto-feed system. The combo kits include Infant Flow LP nCPAP Generator kits and single-limb, heated breathing circuits.

Part #	Description
7772011	Infant Flow LP nCPAP Generator/Circuit F&P MR 730 & 850 Humidifier, 10 per box
7772021	Infant Flow nCPAP Generator/Circuit F&P MR 730 & 850 Humidifier, 20 per box



Infant Flow™ LP nCPAP Starter Kit

The Infant Flow LP nCPAP Starter kit enables your customer to evaluate the headgear fixation. The kit includes:

- 5 Infant Flow LP Generator kits, which include the Infant Flow LP nCPAP generator, a sizing guide, and one each of small, medium and large nasal prongs
- 2 extra large nasal prongs
- 2 extra small nasal prongs
- 2 of each mask sizes (extra small, small, medium, large, extra large)
- 2 of each headgear sizes (extra small, small, small-medium, medium, large and extra large)


Part #	Description
7772010	Infant Flow LP nCPAP Starter Kit



Infant Flow™ LP Headgear

The unique design of the headgear provides a custom fit and secures the generator to headgear. The intuitive “comfort wrap” allows a single clinician to quickly apply the headgear. The headgear is available in six sizes to fit head circumferences from 17 to 42 centimeters. Individually packaged 10 per case.

Part #	Description
777040XS	Infant Flow LP Headgear, Extra Small
777040S	Infant Flow LP Headgear, Small
777040SM	Infant Flow LP Headgear, Small-Medium
777040M	Infant Flow LP Headgear, Medium
777040L	Infant Flow LP Headgear, Large
777040XL	Infant Flow LP Headgear, Extra Large



Infant Flow™ LP Bonnet

The traditional soft cotton bonnets are available in 10 sizes to fit head measurements from 18 cm to 38 cm. The attachment pad in the center of the bonnet safely secures the generator to the bonnet. A color-coded edge enables easy size identification. Bonnets are packaged 10 per case.



Part #	Description
777010	LP Bonnet, Size 000
777012	LP Bonnet, Size 00
777014	LP Bonnet, Size 0
777016	LP Bonnet, Size 1
777018	LP Bonnet, Size 2

Part #	Description
777020	LP Bonnet, Size 3
777022	LP Bonnet, Size 4
777024	LP Bonnet, Size 5
777026	LP Bonnet, Size 6
777028	LP Bonnet, Size 7

Infant Flow™ Generator Straps

The gray straps attach the Infant Flow LP generator to the bonnet. The straps are used to replace straps if they lose their shape. The straps are available 20 per package.



Part #	Description
775728	Infant Flow Generator Straps



GLOBAL
HEADQUARTERS

Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Vyairē Medical Pty Ltd
Level 5, 7 Eden Park Drive
Macquarie Park, NSW, 2113
Australia

vyairē.com

For Latin America, Canada, EU, Switzerland, South Africa, Australia, Russia and Asia Pacific distribution only.

Trademarks are the property of their respective owners.

© 2019 Vyairē Medical, Inc. or one of its affiliates. All rights reserved. Vyairē, the Vyairē Medical logo and Infant Flow LP nCPAP System are trademarks or registered trademarks of Vyairē Medical, Inc. or one of its affiliates. Medical devices class I and IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-INT-1900251