



SalterSTAT Capnography

The signal you need—anytime, anywhere

Striving
to be
**#1 IN PATIENT
SATISFACTION**




Innovative Options for Everywhere Capnography Is Indicated

		Green Series for Luer-Style Connector Applications	Orange Series for Reflective Connector Applications
			
	Cannulas	Sampling Only	✓
		Nasal with O ₂ Delivery	✓
		Oral-Nasal with O ₂ Delivery	✓
	Sampling Lines	Intubated and Inline	✓
Blue Series for Mask Applications	Masks	Oral-Nasal for Higher Flows	✓
		Oral-Nasal for Upper-GI Procedures	✓

Innovative Designs Are Effective at CO₂ Detection and O₂ Delivery


Divided Designs¹

EtCO₂ sampling and O₂ delivery effectiveness at 4 LPM
Fresh Gas Flow



CapnoVue Masks²

Face masks with CapnoVue® technology
provide accurate signals with high O₂ flow rates



Designed to Perform in Reflective Connector Systems




40% less
dead space in facepiece²

SalterSTAT Orange
Oral-Nasal Sampling Line
with O₂ Delivery



Competitive Oral-Nasal
Sampling Line with
O₂ Delivery



Orange Series
for Reflective
Connector
Systems


Designed for Reliability in Reflective Connector Systems



Compatible connector




Hydrophobic, 0.2 micron filter



**Nafion® tubing
for long-term applications**

Orange Series
for Reflective
Connector
Systems

Divided Cannulas for True EtCO₂ Sampling and Effective O₂ at Higher LPM Rates




■ Patented Salter Eyes provide safety apertures to reduce the possibility of occlusions

■ Divided design separates O₂ delivery from CO₂ monitoring

Green Series
for Luer-Style
Connector
Systems

CapnoVue Masks Provide Higher O₂ Flow Rates with Oral-Nasal Sampling



■ Available in standard design and with access port for scopes

■ Flow dynamics separate gases by pulling CO₂ toward vents and sampling lines and away from O₂ port

Blue Series
for Luer-Style
Connector
Systems


SalterSTAT Capnography

The signal you need—anytime, anywhere


Striving
to be
**#1 IN PATIENT
SATISFACTION**



**Orange Series
for Reflective
Connector Systems**



**Blue Series for
Mask Applications**



**Green Series
for Luer-Style
Connector Systems**

Improve intubation safety with Salter's innovative airway management products



**air-Q Masked
Laryngeal Airway**



**IntuBrite
Laryngoscopes**



**Flex-Tip
Endotracheal Tubes**

1. Ebert TJ, Novaljija J, Uhrich TD, Barney JA. The effectiveness of oxygen delivery and reliability of carbon dioxide waveforms: a crossover comparison of 4 nasal cannulae. *Anesth Analg*. 2015;120(2):342-438.

2. Data on file. Salter Labs.

3. Salter Labs Technical Report 006

30 Spur Drive, El Paso, TX 79906
Phone: 1-800-421-0024
International: 1-303-799-1990
Fax: 800-628-4690
E-mail: customercare@salterlabs.com
www.salterlabs.com

SLML-205 Rev A July 2018

SALTER LABS[®]
DESIGN INNOVATIONS FOR
RESPIRATORY AND AIRWAY CARE