



**TRUDELL**  
HEALTHCARE SOLUTIONS

**NISSHA**  
MEDICAL TECHNOLOGIES

HEALTHCARE  
SOLUTIONS

**VERMED**

# SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

## FEATURES



### MINIMAL CLEANUP

Conductive adhesive solid gel requires minimal cleanup.



### COMFORTABLE AND FUNCTIONAL

Pliable substrate conforms to the body.



### CONVENIENT PACKAGING

Patented traction liner prevents tabs from sliding off the liner and ensures packs do not brick together for ease of removal and application.



### HYPOALLERGENIC

Tab electrodes are not made with natural rubber latex.

Vermed SilveRest™ Tab Electrodes for resting EKG provide excellent adhesion and quality traces.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

### INNER BOX

Single use packs are packaged in inner boxes within the case

Latex-free Solid Gel



# SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

Vermed SilveRest™ Tab Electrodes for resting EKG provide excellent adhesion and quality traces. All Vermed tab electrodes meet or exceed the AAMI standard for disposable pre-gelled EKG electrodes.

## A10022 - .984" x .90" / 24.99mm x 22.86mm - Ag/AgCL



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10022	10	100	1000	4000		24 Hours

## A10201 - 1.34" x .90" / 34.036mm x 22.86mm - C/AgCL



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10201	10	100	1000	4000		24 Hours
A10201-500	10	100	500	3000		

## A10023 - 1.125" x 1.438" / 28.575mm x 36.525mm - Ag/AgCL



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10023	10	100	1000	4000		24 Hours
A10023-1000BX	10	100	n/a	1000		
A10023-10	10	10		4000		

## A10202-10 - 1.297" x .864" / 34mm x 23mm - C/AgCL



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10202-10	10	10	500	4000		24 Hours

## A10027 - 0.5" x 0.9 / 12.7mm x 22.86mm - AG/AGCL TinyTab



TinyTab electrodes are specially designed to fit on smaller pediatric patients

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10027	10	100	n/a	1000		24 Hours