

**A10079-SRT** (3 Lead)  
**A10077-SRT** (5 Lead)  
Cloth | Solid Gel

**A10076**  
Hydrocolloid | Foil

**A10080-SRT**  
Cloth | Solid Gel

## FEATURES

### 5 Day WEAR TIME GUARANTEED

**A10079-SRT** (3 Lead) | Cloth | Solid Gel  
**A10077-SRT** (5 Lead) | Cloth | Solid Gel

- Long-term high humidity ECG Electrode
- Hydrocolloid border provides long term attachment

**A10080-SRT** | Cloth | Solid Gel

- Micro size can be used on limbs or on body
- Safely secures to limbs without a strap that can reduce circulation
- Gentle hydrocolloid and hydrogel adhesives

**A10076** | Hydrocolloid | Foil

- Hydrocolloid is gentle on delicate skin
- Provides long term attachment with no skin irritation
- Resists breakdown from moisture or heat
- Leaves no residue on skin

**1.800.669.6905 | Vermed.com**

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

MR Conditional



Latex-free



Solid Gel



# TenderTrode Plus™


## Pre-Wired Electrodes & Thermal Probe Covers

Vermed's TenderTrode Plus™ line of hydrocolloid products are designed for the unique needs of the neonatal market. Each product uses a hydrocolloid instead of aggressive tape adhesives to prevent irritation to tender skin while providing long term attachment. This, coupled with the use of soft conforming materials, ensures patient comfort and product effectiveness in the NICU.


### A10079-SRT - 1" / 25.4MM - ROUND PREWIRED CLOTH ELECTRODE - 3 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10055-SRT	3	3	n/a	60		

### A10077-SRT - 1" / 25.4MM - ROUND PRE-WIRED CLOTH ELECTRODE - 5 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10077-SRT	5	5	n/a	100		

### A10080-SRT - 1" X 0.375" / 25.4MM X 9.525MM - RECT. PRE-WIRED CLOTH ELECTRODE - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10058-SRT	3	3	n/a	60		

### A10076 - 1" / 25.4MM - HYDROCOLLOID PROBE COVER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10076	1	50	n/a	500	

## ADDITIONAL FEATURES

- Save nursing time and resources with fewer electrode changes
- Safe on neonatal skin
- Designed specifically for NICU
- Radiolucent Lead

## Clinical Feedback

Recent clinical feedback indicates a **98% reduction in faulty incubator and monitor alarms** tripped when standard electrodes peel off in high heat and humidity.

TenderTrode™ Plus is now **leading the way** in neonatal patient care. **Significant cost savings** are evident due to limited replacement, a function of our patented hydrocolloid adhesive border. These electrodes have a 5 day wear time - Guaranteed.

Reducing the removal and replacement of standard monitoring electrodes **ensures the body is focused on essential growth and repair.** Reduced monitor alarms ensures the comfort of patients and parents alike while reducing nursing staff workload.