


Pre-filled 0.9% Sodium Chloride I.V. Flush Syringes



TRUDELL
HEALTHCARE SOLUTIONS

TRU TEST
A Tradition of Trust.
Une tradition de confiance.

The **TruTest® 0.9% Sodium Chloride I.V. Flush Syringes** have an optimized syringe assembly to improve patient outcomes.

- Needle free design eliminates risk of an inadvertent needlestick
- Single use design to help reduce risk of medication errors and contamination
- Distinct and prominent syringe labelling allows for easy identification of syringe size and volume

Key Features:

- Luer lock tip
- Latex and DEHP free
- Disposable
- 12 mL syringe offered in **3 mL**, **5 mL** and **10 mL** volumes



10 mL Pre-filled I.V. Flush Syringe



5 mL Pre-filled I.V. Flush Syringe



3 mL Pre-filled I.V. Flush Syringe

Ordering Information:

Item Number	Description	Unit of Measure
T161203-BP	Syringe, I.V. Flush, 0.9% Sodium Chloride, 3mL FILL in 12mL Syringe	Case/400
T161205-BP	Syringe, I.V. Flush, 0.9% Sodium Chloride, 5mL FILL in 12mL Syringe	Case/400
T161210-BP	Syringe, I.V. Flush, 0.9% Sodium Chloride, 10mL FILL in 12mL Syringe	Case/400



HEALTHPRO

Contract Number:

PH401RAF_PHCL

For more information, please contact your local THS Territory Manager!

1-800-265-5494

customerservice@trudellhs.com

www.trudellhs.com