

Precision Medical offers a variety of vacuum regulator solutions to meet the needs of your customer. Please reference the chart on page 2 for the full list of product numbers.

## Most Common Products:

Product Code:	Description:
PM3303	Cont-Int, DISS H.T., Tubing Nipple, 0-200 mmHg
PM3314	Cont-Int, DISS H.T., DISS Male, 0-200 mmHg
PM3103	Cont, DISS H.T., Tubing Nipple, 0-200 mmHg
PM3114	Cont, DISS H.T., DISS Male, 0-200 mmHg

## Most Common Accessory: VacTrap

Product Code:	Description:
PM9000	Base Model
PM9001	Ohmeda Style Trap Bottle Connection
PM9002	DISS Hex
PM9003	DISS Hand Tight

- Shatter Resistant Trogamid Bottle
- Lightweight Nylon Float
- Positive Shut-Off Disc Seal
- Designed For Easy Disassembly



## There are 4 important questions to ask your customer:

### 1. Continuous or Continuous-Intermittent model?

#### Continuous & Continuous-Intermittent Models:

Continuous	Continuous-Intermittent
PM3000 PM3100	PM3300
PM3500 PM3600	PM3400




### 2. How much suction is required? (Gauge Range in mmHg)

#### Continuous Models

	Gauge Range:
PM3000 & PM3100	0-200 mmHg and "Full Vacuum"
PM3000HV & PM3100HV	0-300 mmHg and "Max Vacuum"
PM3500	0-150 mmHg
PM3600	0-760 mmHg

#### Continuous-Intermittent Models


	Gauge Range:
PM3300	0-200 mmHg and "Full Vacuum"
PM3300HV	0-300 mmHg and "Max Vacuum"
PM3400	0-150 mmHg



### 3. What is the fitting at the back of the regulator? (Vacuum Source)

### 4. What is the fitting at the bottom of the regulator? (Inlet)

		VACUUM SOURCE CONNECTION (on back of unit)										
		1/8" N.P.T.F.	DISS HEX	DISS H.T.	OHMEDA	CHEMETRON	OXEQUIP (Old Style)	OXEQUIP (Medstar)	PURITAN-BENNETT	SCHRADER	TUBING NIPPLE	90° TUBING NIPPLE
BOTTOM CONNECTION (inlet)		PM3001	PM3002	PM3003	PM3004	PM3005	PM3006	PM3007	PM3008	PM3009	PM3010	PM3011
		PM3101	PM3102	PM3103	PM3104	PM3105	PM3106	PM3107	PM3108	PM3109	PM3110	PM3111
		PM3501	PM3502	PM3503	PM3504	PM3505	PM3506	PM3507	PM3508	PM3509	PM3510	PM3511
		PM3601	PM3602	PM3603	PM3604	PM3605	PM3606	PM3607	PM3608	PM3609	PM3610	PM3611
		PM3301	PM3302	PM3303	PM3304	PM3305	PM3306	PM3307	PM3308	PM3309	PM3310	PM3311
		PM3401	PM3402	PM3403	PM3404	PM3405	PM3406	PM3407	PM3408	PM3409	PM3410	PM3411
		PM3012	PM3013	PM3014	PM3015	PM3016	PM3017	PM3018	PM3019	PM3020	PM3021	PM3022
		PM3112	PM3113	PM3114	PM3115	PM3116	PM3117	PM3118	PM3119	PM3120	PM3121	PM3122
		PM3512	PM3513	PM3514	PM3515	PM3516	PM3517	PM3518	PM3519	PM3520	PM3521	PM3522
		PM3612	PM3613	PM3614	PM3615	PM3616	PM3617	PM3618	PM3619	PM3620	PM3621	PM3622
		PM3312	PM3313	PM3314	PM3315	PM3316	PM3317	PM3318	PM3319	PM3320	PM3321	PM3322
		PM3412	PM3413	PM3414	PM3415	PM3416	PM3417	PM3418	PM3419	PM3420	PM3421	PM3422
		PM3023	PM3024	PM3025	PM3026	PM3027	PM3028	PM3029	PM3030	PM3031	PM3032	PM3033
		PM3123	PM3124	PM3125	PM3126	PM3127	PM3128	PM3129	PM3130	PM3131	PM3132	PM3133
		PM3523	PM3524	PM3525	PM3526	PM3527	PM3528	PM3529	PM3530	PM3531	PM3532	PM3533
		PM3623	PM3624	PM3625	PM3626	PM3627	PM3628	PM3629	PM3630	PM3631	PM3632	PM3633
		PM3323	PM3324	PM3325	PM3326	PM3327	PM3328	PM3329	PM3330	PM3331	PM3332	PM3333
		PM3423	PM3424	PM3425	PM3426	PM3427	PM3428	PM3429	PM3430	PM3431	PM3432	PM3433
		PM3034	PM3035	PM3036	PM3037	PM3038	PM3039	PM3040	PM3041	PM3042	PM3043	PM3044
		PM3134	PM3135	PM3136	PM3137	PM3138	PM3139	PM3140	PM3141	PM3142	PM3143	PM3144
PM3534		PM3535	PM3536	PM3537	PM3538	PM3539	PM3540	PM3541	PM3542	PM3543	PM3544	
PM3634		PM3635	PM3636	PM3637	PM3638	PM3639	PM3640	PM3641	PM3642	PM3643	PM3644	
	PM3334	PM3335	PM3336	PM3337	PM3338	PM3339	PM3340	PM3341	PM3342	PM3343	PM3344	
	PM3434	PM3435	PM3436	PM3437	PM3438	PM3439	PM3440	PM3441	PM3442	PM3443	PM3444	
	PM3000	PM3046	PM3047	PM3048	PM3049	PM3050	PM3051	PM3052	PM3053	PM3054	PM3055	
	PM3100	PM3146	PM3147	PM3148	PM3149	PM3150	PM3151	PM3152	PM3153	PM3154	PM3155	
	PM3500	PM3546	PM3547	PM3548	PM3549	PM3550	PM3551	PM3552	PM3553	PM3554	PM3555	
	PM3600	PM3646	PM3647	PM3648	PM3649	PM3650	PM3651	PM3652	PM3653	PM3654	PM3655	
	PM3300	PM3346	PM3347	PM3348	PM3349	PM3350	PM3351	PM3352	PM3353	PM3354	PM3355	
	PM3400	PM3446	PM3447	PM3448	PM3449	PM3450	PM3451	PM3452	PM3453	PM3454	PM3455	



**TRUCELL**  
HEALTHCARE SOLUTIONS